

# Measure breath nitric oxide for airway inflammation with the **NObreath®** FeNO Device



*Aids in the diagnosis and management of asthma, one breath at a time.*

[www.bedfont.com](http://www.bedfont.com)





**For more information on this product call us now on**

## Contents

Fractional Exhaled Nitric Oxide (FeNO)	4
Benefits of performing FeNO tests	4
NObreath® features	6
Measuring FeNO with NObreath®	7
Consumables	7
Technical specification	8
FeNOchart™	9
NObreath® Forum	9
Interpretation chart	10-11
References	11

## Fractional Exhaled Nitric Oxide (FeNO)

Fractional exhaled Nitric Oxide (FeNO) is a good marker for eosinophilic airway inflammation, and is considered to be a good indicator of corticosteroid response<sup>1</sup>.

The production of nitric oxide is often found to be higher in inflammatory conditions such as asthma; therefore FeNO monitoring can be used for the detection and management of such conditions<sup>2</sup>, but also to differentiate between COPD, ACOS and other interstitial lung diseases that are not assessed by other means, such as lung function<sup>3</sup>.

Nitric oxide measurement is not intended as a stand-alone method for diagnosis and should be used in conjunction with other evaluation methods and tests<sup>4</sup>.

FeNO measurement is a simple, rapid, highly reproducible, and non-invasive method of airway inflammation assessment, which until now, has been an expensive test to deliver in everyday practice<sup>5</sup>.


### Benefits of performing FeNO tests:

- Non-invasive, quick and easy to perform<sup>5</sup>.
- Aids in asthma management, assisting with the correct prescription and making monitored adjustments.
- Shows patient adherence to treatment<sup>6</sup>.
- Aids in identifying good and poor adherence to corticosteroid treatment<sup>1</sup>.
- Good indicator of corticosteroid response<sup>1</sup>.
- Shown to be superior to the majority of conventional tests of lung function, such as peak flow recording and spirometry<sup>5</sup>.
- Aids in differentiating between allergic (eosinophilic) and non-allergic asthma<sup>7</sup>.



# NObreath® features

An ergonomic design, fully-portable and incorporated with antimicrobial technology for optimum infection control.




\*Subject to correct use, maintenance and service. Tested up to 29,000 tests.

For more information on this product call us now on

# Measuring FeNO with NObreath®

IT'S AS EASY AS:



## Consumables

### NObreath® Mouthpiece

The NObreath® uses a single-patient use mouthpiece, which contains an integrated infection control filter that removes and traps > 99% of airborne bacteria and > 98% of viruses<sup>8</sup>.

Dimensions	Approx. 180 mm x 28 mm x 22 mm
Weight	Approx. 14 g
Material	Polypropylene
Shelf life	5 Years



## Technical specification

Concentration range		5 - 500 ppb
Display		Full colour touchscreen
Detection principle		Electrochemical sensor
Repeatability		± 5 ppb of measured value ≤ 50 ppb ± 10% of measured value > 50 ppb
Accuracy		± 5 ppb of measured value ≤ 50 ppb ± 10% of measured value > 50 ppb
Power	NObreath® Device	1 x main rechargeable Li-ion battery - Approx. 100 uses on fully charged battery 2 x Li-ion coin cell battery- Approx. 5 years Input: 5 V, 0.5 A
	NObreath® Dock	Mains powered Input: 5 V, 0.5 A Output: 5 V, 0.5 A
	Plug	Input: 100 - 240 V ~ 50/60 Hz., 0.2 A Output: 5.0 V, 1.0 A
T <sub>90</sub> response time		≤ 10 seconds
Temperature	Operating	15 - 30°C
	Storage/transport	0 - 50°C
Humidity	Operating	20 - 80% RH (non-condensing)
	Storage/transport	5 - 95% RH (non-condensing)
Operating/storage/transport pressure		800 – 1080 mbar
Sensor operating life		5 years (Subject to servicing)
Sensor sensitivity		1 ppb
Sensor drift		< 5% per annum
Dimensions		Approx. 90 mm x 159 mm x 59 mm
Weight		Approx. 400 g
Materials	NObreath® Device	Case: polycarbonate/ABS blend with antimicrobial technology
	NObreath® Dock	
Breath test time		Adult: 12 seconds Child: 10 seconds Ambient: 30 seconds
Warm-up time		≤ 60 seconds
Maximum ambient operating level		350 ppb NO
CO cross interference		45 ppm ≤ 17.6 ppb

**NOTE:** Exhaled flow during FeNO measurement at 50 ml / sec ± 10% at 10 cm H<sub>2</sub>O

## FeNOchart™

FeNOchart™ is a free patient management software, available with every NObreath®. FeNOchart™ enables you to track patients' progress, view live readings, download results plus much more.



**FREE FeNOchart™ patient management software.**

## NObreath® Forum

Purchasing a NObreath® entitles you to free membership of the NObreath® forum. The NObreath® forum is an international, invitation-only platform where professionals using the Bedfont® NObreath® FeNO device can communicate, share experiences and knowledge, and ask for other professional opinions. There is no cost or obligation to participate and membership is free when you purchase a NObreath®.



# Using FeNO to assist diagnosis

Measuring airway inflammation with NObreath® can help monitor the effectiveness

## Aid in diagnosis using the NObreath® FeNO device

FeNO (ppb) Levels	LOW < 25 ppb (< 20 ppb in children)	INTERMEDIATE 25 - 50 ppb (20 - 35 ppb in children)	HIGH > 50 ppb (> 35 ppb in children) or rise in FeNO of > 40% from previously stable levels
<b>Symptomatic</b> (chronic cough and/or wheeze and/or shortness of breath during past 6 wk)	**Eosinophilic airway inflammation unlikely  Alternative diagnosis  Unlikely to benefit from ICS	Be cautious  Evaluate clinical context  Monitor change in FeNO over time	**Eosinophilic airway inflammation present  Likely to benefit from ICS

## Alternative considerations (if Allergic Asthma has been dismissed)<sup>2</sup>

- Non-Allergic Asthma
- Chronic cough
- Vocal Chord Dysfunction
- GERD

Nitric oxide measurement is not intended as a stand-alone method for diagnosis and should be used in conjunction with other evaluation methods and tests<sup>4</sup>.



# & management of Asthma

of medication and can be used to predict the risk of Asthma attacks\*.

## Monitoring (in patients with diagnosed asthma) using the NObreath® FeNO device

FeNO (ppb) Levels	LOW < 25 ppb (< 20 ppb in children)	INTERMEDIATE 25 - 50 ppb (20 - 35 ppb in children)	HIGH > 50 ppb (> 35 ppb in children) or rise in FeNO of > 40% from previously stable levels
<b>Symptomatic</b> (chronic cough and/or wheeze and/or shortness of breath during past 6 wk)	Possible alternative diagnosis  Unlikely to benefit from increase in ICS	Persistent allergen exposure  Inadequate ICS dose  Poor adherence  Steroid resistance	Persistent allergen exposure  Poor adherence or inhaler technique  Inadequate ICS dose  Risk of Exacerbation  Steroid resistance
<b>Symptoms Absent</b>	Adequate ICS dose  Good adherence  ICS taper	Adequate ICS dosing  Good adherence  Monitor Change in FeNO	ICS withdrawal or dose reduction may result in relapse  Poor adherence or inhaler technique

## References

1. Price D, Ryan D, Burden A, Von Ziegenweidt J, Gould S, Freeman D et al. Using fractional exhaled nitric oxide (FeNO) to diagnose steroid-responsive disease and guide asthma management in routine care. *Clinical and Translational Allergy*. 2013;3(1).
2. Saito J, Gibeon D, Macedo P, Menzies-Gow A, Bhavsar P, Chung K. Domiciliary diurnal variation of exhaled nitric oxide fraction for asthma control. 2017.
3. ATS/ERS Recommendations for Standardized Procedures for the Online and Offline Measurement of Exhaled Lower Respiratory Nitric Oxide and Nasal Nitric Oxide, 2005; *American Journal of Respiratory and Critical Care Medicine*; vol. 171: 912-930;2005.
4. Correlation of Exhaled Nitric Oxide, Spirometry and Asthma Symptoms: *Journal of Asthma*: Vol 42, No 10 [Internet]. Tandfonline.com. 2017 [cited 15 March 2017]. Available from: <http://www.tandfonline.com/doi/abs/10.1080/02770900500371344>.
5. Chen F, Liao H, Huang X, Xie C. Importance of fractional exhaled nitric oxide in diagnosis of bronchiectasis accompanied with bronchial asthma. *Journal of Thoracic Disease*. 2016;8(5):992-999.
6. Beck-Ripp J, Griese M, Arenz S, Koring C, Pasqualoni B, Buefler P. Changes of exhaled nitric oxide during steroid treatment of childhood asthma. *Eur Respir J* 2002;19:1015–1019.
7. Coumou HBel E. Improving the diagnosis of eosinophilic asthma [Internet]. Taylor and Francis online. 2017 [cited 15 March 2017]. Available from: <http://www.tandfonline.com/doi/full/10.1080/17476348.2017.1236688>.
8. Public Health England. An Evaluation of Filtration Efficiencies Against Bacterial and Viral Aerosol Challenges. Salisbury: Public Health England; 2020.

\* FeNO is not a definitive indication of asthma and should be used in conjunction with (but not limited to) spirometry, patient history, symptoms.

\*\* Allergic = Eosinophilic / Non- Allergic = Non-Eosinophilic.



Contact Bedfont® or one of our worldwide **NObreath®** distributors for a free demonstration.

[www.bedfont.com](http://www.bedfont.com)  
Tel: +44 (0)1622 851122  
Email: [ask@bedfont.com](mailto:ask@bedfont.com)

***World leaders in breath analysis.***

Visit [www.bedfont.com/resources](http://www.bedfont.com/resources) to view this document in other languages.



**Bedfont® Scientific Ltd.**  
Station Road, Harrietsham, Maidstone,  
Kent, ME17 1JA England  
Tel: +44 (0)1622 851122 Fax: +44 (0)1622 854860  
Email: [ask@bedfont.com](mailto:ask@bedfont.com) Web: [www.bedfont.com](http://www.bedfont.com)



**Emergo Europe B.V.**  
Westervoortsedijk 60  
6827 AT Arnhem  
The Netherlands.

© Bedfont® Scientific Limited 2024

Issue 11 - October 2024, Part No: MKT503  
Bedfont® Scientific Limited reserves the right to change or update this literature without prior notice.  
Registered in: England and Wales. Registered No: 1289798



MD 502905